

# Tidewater Regional IV Box Layout

|                                                                                                                                                              |  |                                                                                                                                                                                              |                                   |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------|--|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------|
| 20 Gauge Needles 1" (3) Safety<br>22 Gauge Needles 1.5" (3) Safety<br>25 Gauge Needles 5/8" (3) Safety                                                       |  | Medication Labels (x4)<br><br>T.B. Syringe (x3)<br><br>Blunt Needles 18g (x3) Safety<br><br>Bioclusive Dressings (x3)<br><br>Filter Straw<br><br>10 ml Syringe (x3)<br><br>3 ml Syringe (x3) | Sharps Container Alcohol Preps(8) |
| 10 ml Sterile Water (x2)<br>10ml Sterile Saline 0.9% (x2) Pre-fill<br>Roll 1" Tape (x2)<br>Tourniquets(x2)<br>Nitrostat Tabs<br>Albuterol (x3) Atrovent (x1) |  |                                                                                                                                                                                              |                                   |
| Injectable Ports(4) Aspirin (8) Chewable 81mg ea. Band-aids                                                                                                  |  |                                                                                                                                                                                              |                                   |

Top Shelf

|                                                                                                                                                                         |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Mini-drip IV admin set (60 gtts/cc)<br>24 Gauge Angicaths (x3) 3/4" Safety                                                                                              |
| 22 Gauge Angicaths (x3) 1" Safety                                                                                                                                       |
| 20 Gauge Angicaths (x3) 1" Safety<br>Intraosseous Needle (15/16 ga.)                                                                                                    |
| 18 Gauge Angicaths (x3) 1 1/2" Safety                                                                                                                                   |
| Macro-drip IV admin set<br>(10 or 15 gtts/ml)<br>16 Gauge Angicaths (x3) 1 1/2" Safety<br>14 Gauge Angicaths (x2) 2" Non-Safety<br>14g shall not be "protectocath" type |

Middle Shelf

|                             |                      |
|-----------------------------|----------------------|
| Hand-held Nebulizer         | 4-way Stopcocks (x4) |
| Macro-drip IV admin set(x2) | 100 mL 0.9% NS (x1)  |
| (10-15 gtts/mL)             | 1000 mL 0.9% NS (x3) |
|                             |                      |
| Mini-drip IV admin set(x2)  |                      |
| (60 gtts/mL)                |                      |
| 50 mL Syringe (2)           | 9" Arm Boards (x2)   |
|                             | 250 mL 0.9% NS (x1)  |

Bottom of Box



Janaury 2010